

## RAMP1 rabbit pAb

Cat No.: ES11932

For research use only

## Overview

Product Name RAMP1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

**Specificity** RAMP1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Receptor activity-modifying protein 1

(Calcitonin-receptor-like receptor activity-modifying

protein 1) (CRLR activity-modifying protein 1)

Gene Name RAMP1

**Cellular localization** Membrane; Single-pass type I membrane protein. **Purification** The antibody was affinity-purified from rabbit

ification The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band16kDHuman Gene ID10267

**Human Swiss-Prot Number** 

**Alternative Names** 

**Background** The protein encoded by this gene is a member of

060894

the RAMP family of single-transmembrane-domain proteins, called receptor (calcitonin) activity modifying proteins (RAMPs). RAMPs are type I transmembrane proteins with an extracellular N terminus and a cytoplasmic C terminus. RAMPs are

required to transport calcitonin-receptor-like



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com



receptor (CRLR) to the plasma membrane. CRLR, a receptor with seven transmembrane domains, can function as either a calcitonin-gene-related peptide (CGRP) receptor or an adrenomedullin receptor, depending on which members of the RAMP family are expressed. In the presence of this (RAMP1) protein, CRLR functions as a CGRP receptor. The RAMP1 protein is involved in the terminal glycosylation, maturation, and presentation of the CGRP receptor to the cell surface. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided b

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

