

ZRAN1 rabbit pAb

Cat No.: ES12076

For research use only

Overview

Product Name ZRAN1 rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human ZRAN1 AA

range: 255-305

Specificity This antibody detects endogenous levels of ZRAN1

at Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name ZRAN1

Gene Name ZRANB1 TRABID

Cellular localization Cytoplasm . Nucleus . Enriched in punctate

localization in the cytoplasm. .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 54764 Human Swiss-Prot Number Q9UGI0


Alternative Names

+86-27-59760950

Background







Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night

