

# ZN879 rabbit pAb

Cat No.:ES12094

For research use only


## Overview

|                          |                                                                                                                       |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name             | ZN879 rabbit pAb                                                                                                      |
| Host species             | Rabbit                                                                                                                |
| Applications             | WB                                                                                                                    |
| Species Cross-Reactivity | Human; Mouse                                                                                                          |
| Recommended dilutions    | WB 1:500-2000                                                                                                         |
| Immunogen                | Synthesized peptide derived from human ZN879 AA range: 446-496                                                        |
| Specificity              | This antibody detects endogenous levels of ZN879 at Human/Mouse                                                       |
| Formulation              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage                  | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                                    |
| Protein Name             | ZN879                                                                                                                 |
| Gene Name                | ZNF879                                                                                                                |
| Cellular localization    | Nucleus .                                                                                                             |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality                | Polyclonal                                                                                                            |
| Concentration            | 1 mg/ml                                                                                                               |
| Observed band            | 60kD                                                                                                                  |
| Human Gene ID            | 345462                                                                                                                |
| Human Swiss-Prot Number  | B4DU55                                                                                                                |
| Alternative Names        | Zinc finger protein 879                                                                                               |
| Background               |                                                                                                                       |





**ELK Biotechnology**



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.