

ZN358 rabbit pAb

Cat No.: ES12169

For research use only

Overview

Product Name ZN358 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1:500-2000

Immunogen Synthesized peptide derived from human ZN358 AA

range: 436-486

Specificity This antibody detects endogenous levels of ZN358 at

Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name ZN358
Gene Name ZNF358
Cellular localization Nucleus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 60kD
Human Gene ID 140467
Human Swiss-Prot Number Q9NW07


Alternative Names Zinc finger protein 358

Background

+86-27-59760950







Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

