

ZN350 rabbit pAb

Cat No.: ES12170

For research use only

Overview

Product Name ZN350 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human ZN350 AA

range: 346-396

Specificity This antibody detects endogenous levels of ZN350 at

Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name ZN350

Gene Name ZNF350 ZBRK1

Cellular localization Nucleus. Nucleus matrix. Associated with the

nuclear matrix.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

Observed band


Human Gene ID 59348 Human Swiss-Prot Number Q9GZX5

Alternative Names

Background







Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

