

ZN140 rabbit pAb

Cat No.: ES12192

For research use only

Overview

Product Name ZN140 rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human ZN140 AA

range: 222-272

Specificity This antibody detects endogenous levels of ZN140 at

Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name ZN140
Gene Name ZNF140
Cellular localization Nucleus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 7699 Human Swiss-Prot Number P52738


Alternative Names

+86-27-59760950

Background







Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night

