

ZFP57 rabbit pAb

Cat No.:ES12208

For research use only

Overview

Product Name	ZFP57 rabbit pAb	
Host species	Rabbit	
Applications	WB	
Species Cross-Reactivity	Human; Mouse;Rat	
Recommended dilutions	WB 1: 500-2000	
Immunogen	Synthesized peptide derived from human ZFP57 AA range: 373-423	
Specificity	This antibody detects endogenous levels of ZFP57 at Human/Mouse/Rat	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	ZFP57	
Gene Name	ZFP57 C6orf40 ZNF698	
Cellular localization	Nucleus . Binds various differentially methylated regions (DMR)	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration Observed band	1 mg/ml	
Human Gene ID	346171	
Human Swiss-Prot Number Alternative Names	Q9NU63	
Background	The protein encoded by this gene is a zinc finger protein containing a KRAB domain. Studies in mouse suggest that this protein may function as a transcriptional repressor. Mutations in this gene have been associated with transient neonatal diabetes mellitus type 1 (TNDM1).[provided by RefSeq, Sep 2009],	




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C