

ZDHC6 rabbit pAb

Cat No.:ES12218

For research use only

Overview

Product Name	ZDHC6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human ZDHC6 AA range: 160-210
Specificity	This antibody detects endogenous levels of ZDHC6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ZDHC6
Gene Name	ZDHHC6 ZNF376
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass
	membrane protein . When not palmitoylated,
	accumulates to dot-like structures in the
	endoplasmic reticulum (PubMed:28826475)
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	64429
Human Swiss-Prot Number	Q9H6R6
Alternative Names	
Background	




ELKbio@ELKbiotech.com

www.elkbiotech.com

www.eikbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Western blot analysis of lysates from HuvEc cells, primary antibody was diluted at 1:1000, 4°over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C