

XYLK rabbit pAb

Cat No.: ES12273

For research use only

Overview

Product Name XYLK rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human XYLK AA

range: 121-171

Specificity This antibody detects endogenous levels of XYLK at

Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name XYLK

Gene Name FAM20B KIAA0475

Cellular localization Golgi apparatus membrane ; Single-pass type II

membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 9917 Human Swiss-Prot Number 075063


Alternative Names

+86-27-59760950

Background







Western blot analysis of lysates from Hela cells, primary antibody was diluted at 1:1000, 4° over night

