

TM230 rabbit pAb

Cat No.:ES12658

For research use only


Overview

Product Name	TM230 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human TM230 AA range: 32-82
Specificity	This antibody detects endogenous levels of TM230 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	TM230
Gene Name	TMEM230 C20orf30 HSPC274 UNQ2432/PRO4992
Cellular localization	Membrane ; Multi-pass membrane protein . Golgi apparatus, trans-Golgi network . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle . Early endosome . Recycling endosome . Late endosome . Cytoplasmic vesicle, autophagosome .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	29058
Human Swiss-Prot Number	Q96A57
Alternative Names	
Background	





ELK Biotechnology



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.