

## TBC3G rabbit pAb

Cat No.: ES12800

For research use only

## Overview

Product Name TBC3G rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human TBC3G AA

range: 231-281

**Specificity** This antibody detects endogenous levels of TBC3G

at Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name TBC3G
Gene Name TBC1D3G

Cell membrane; Lipid-anchor. Associated with lipid

rafts.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 654341 Human Swiss-Prot Number Q6DHY5


**Alternative Names** 

+86-27-59760950

Background







Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

