

## T106B rabbit pAb

Cat No.: ES12849

For research use only

## Overview

Product Name T106B rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human; Mouse; Rat Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human T106B AA

range: 72-122

**Specificity** This antibody detects endogenous levels of T106B at

Human/Mouse/Rat

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name T106B
Gene Name TMEM106B

Cellular localization Late endosome membrane; Single-pass type II

membrane protein. Lysosome membrane;

Single-pass type II membrane protein. Membrane;

Lipid-anchor .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 54664 Human Swiss-Prot Number Q9NUM4


**Alternative Names** 

+86-27-59760950

**Background** 







Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

