

STXB2 rabbit pAb

Cat No.:ES12885

For research use only

Overview

Product Name	STXB2 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA;IHC
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Immunogen	Synthesized peptide derived from human STXB2 AA range: 392-442
Specificity	This antibody detects endogenous levels of STXB2 at
	Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	STXB2
Gene Name	STXBP2 UNC18B
Cellular localization	cytosol, plasma membrane, integral component of
	membrane, apical plasma membrane, SNARE
	complex,specific granule,azurophil granule,zymogen
	granule membrane, cytolytic granule, extracellular
	exosome, tertiary granule,
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	6813
Human Swiss-Prot Number	Q15833
Alternative Names	
Background	This gene encodes a member of the
	STXBP/unc-18/SEC1 family. The encoded protein is
	involved in intracellular trafficking, control of SNARE
	(soluble NSF attachment protein receptor) complex
	assembly, and the release of cytotoxic granules by



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




natural killer cells. Mutations in this gene are associated with familial hemophagocytic lymphohistiocytosis. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Jan 2013],

Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human oophoroma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C