

# SO6A1 rabbit pAb


Cat No.:ES13025

For research use only


## Overview

|                          |                                                                                                                       |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name             | SO6A1 rabbit pAb                                                                                                      |
| Host species             | Rabbit                                                                                                                |
| Applications             | WB;IHC                                                                                                                |
| Species Cross-Reactivity | Human;Rat;Mouse;                                                                                                      |
| Recommended dilutions    | WB 1:500-2000;IHC-p 1:50-300                                                                                          |
| Immunogen                | Synthesized peptide derived from human SO6A1 AA range: 213-263                                                        |
| Specificity              | This antibody detects endogenous levels of SO6A1 at Human                                                             |
| Formulation              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage                  | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                                    |
| Protein Name             | SO6A1                                                                                                                 |
| Gene Name                | SLCO6A1 OATP6A1 SLC21A19                                                                                              |
| Cellular localization    | Cell membrane; Multi-pass membrane protein.                                                                           |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality                | Polyclonal                                                                                                            |
| Concentration            | 1 mg/ml                                                                                                               |
| Observed band            |                                                                                                                       |
| Human Gene ID            | 133482                                                                                                                |
| Human Swiss-Prot Number  | Q86UG4                                                                                                                |
| Alternative Names        |                                                                                                                       |
| Background               |                                                                                                                       |





Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

