



# S35E1 rabbit pAb

Cat No.:ES13237

For research use only


## Overview

|                                 |   |
|---------------------------------|---|
| <b>Product Name</b>             | S35E1 rabbit pAb  |
| <b>Host species</b>             | Rabbit  |
| <b>Applications</b>             | WB  |
| <b>Species Cross-Reactivity</b> | Human; Mouse;Rat  |
| <b>Recommended dilutions</b>    | WB 1: 500-2000  |
| <b>Immunogen</b>                | Synthesized peptide derived from human S35E1 AA range: 257-307  |
| <b>Specificity</b>              | This antibody detects endogenous levels of S35E1 at Human/Mouse/Rat   |
| <b>Formulation</b>              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Storage</b>                  | Store at -20°C. Avoid repeated freeze-thaw cycles.  |
| <b>Protein Name</b>             | S35E1   |
| <b>Gene Name</b>                | SLC35E1 PSEC0038  |
| <b>Cellular localization</b>    | Membrane ; Multi-pass membrane protein .  |
| <b>Purification</b>             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Clonality</b>                | Polyclonal  |
| <b>Concentration</b>            | 1 mg/ml   |
| <b>Observed band</b>            |   |
| <b>Human Gene ID</b>            | 79939   |
| <b>Human Swiss-Prot Number</b>  | Q96K37  |
| <b>Alternative Names</b>        |   |
| <b>Background</b>               |   |





**ELK Biotechnology**



Western blot analysis of lysates from SH-SY5Y cells,  
primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

[ELKbio@ELKbiotech.com](mailto:ELKbio@ELKbiotech.com)

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.