

RRP44 rabbit pAb

Cat No.:ES13331

For research use only

Overview

Product Name	RRP44 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human RRP44 AA
	range: 262-312
Specificity	This antibody detects endogenous levels of RRP44 at
	Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	RRP44
Gene Name	DIS3 KIAA1008 RRP44
Cellular localization	Cytoplasm . Nucleus, nucleolus . Nucleus,
	nucleoplasm . Nucleus . Predominantly located in
	the nucleus (PubMed:20531386). According to
	PubMed:12429849, found in the nucleolus
	(PubMed:12429849). According to
	PubMed:20531386, excluded from nucleolus
	support
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	105kD
Human Gene ID	22894
Human Swiss-Prot Number	Q9Y2L1
Alternative Names	Exosome complex exonuclease RRP44 (EC 3.1.13)
	(EC 3.1.26) (Protein DIS3 homolog) (Ribosomal
	RNA-processing protein 44)
Deelewaynal	

Background




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkl

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C