

RCAN3 rabbit pAb

Cat No.: ES13436

For research use only

Overview

Product Name RCAN3 rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human RCAN3 AA

range: 86-136

Specificity This antibody detects endogenous levels of RCAN3

at Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name RCAN3

Gene Name RCAN3 DSCR1L2

Cellular localization intracellular,cytoplasm,

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 11123 Human Swiss-Prot Number Q9UKA8


Alternative Names

+86-27-59760950

Background







Western blot analysis of lysates from A549 cells, primary antibody was diluted at 1:1000, 4° over night

