

## RASL3 rabbit pAb

Cat No.: ES13458

For research use only

## Overview

Product Name RASL3 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human RASL3 AA

range: 297-347

**Specificity** This antibody detects endogenous levels of RASL3 at

Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name RASL3
Gene Name RASAL3

Cytoplasm . Cytoplasm, cell cortex .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 


Human Gene ID 64926 Human Swiss-Prot Number Q86YV0

**Alternative Names** 

**Background** 







Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night

