

## **QPCT** rabbit pAb

Cat No.: ES13831

For research use only

## Overview

Product Name QPCT rabbit pAb

Host species Rabbit
Applications WB;IHC

Species Cross-Reactivity Human; Mouse

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300

Immunogen Synthesized peptide derived from human QPCT AA

range: 233-283

**Specificity** This antibody detects endogenous levels of QPCT at

Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name QPCT
Gene Name QPCT
Cellular localization Secreted.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 25797 Human Swiss-Prot Number Q16769

**Alternative Names** 

Background This gene encodes human pituitary glutaminyl

cyclase, which is responsible for the presence of pyroglutamyl residues in many neuroendocrine peptides. The amino acid sequence of this enzyme is 86% identical to that of bovine glutaminyl cyclase.


[provided by RefSeq, Jul 2008],




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human cervical carcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

