

PZR (phospho-Tyr241) rabbit pAb

Cat No.: ES13839

For research use only

Overview

Product Name PZR (phospho-Tyr241) rabbit pAb

Host species Rabbit

Applications WB;ELISA;IHC Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000 **Immunogen** Synthesized phosho peptide around human PZR

(Tyr241)

Specificity This antibody detects endogenous levels of

Human Mouse Rat PZR (phospho-Tyr241)

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name PZR (Tyr241)

Gene Name MPZL1 PZR UNQ849/PRO1787

Cellular localizationMembrane ; Single-pass type I membrane protein .PurificationThe antibody was affinity-purified from rabbitPurificationPurification antibody was affinity chromatography using a stress or the purification and the purification

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 30kD
Human Gene ID 9019
Human Swiss-Prot Number 095297


Alternative Names Myelin protein zero-like protein 1 (Protein

zero-related)

Background







Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

