

PPIG (phospho-Ser376) rabbit pAb

Cat No.: ES13995

For research use only

Overview

Product Name PPIG (phospho-Ser376) rabbit pAb

Host species Rabbit
Applications WB;IHC

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300

Immunogen Synthesized phosho peptide around human PPIG

(Ser376)

Specificity This antibody detects endogenous levels of

Human Mouse Rat PPIG (phospho-Ser376)

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name PPIG (Ser376)

Gene Name PPIG

Cellular localization Nucleus matrix . Nucleus speckle . Colocalizes with

RNA splicing factors at nuclear speckles. .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 82kD
Human Gene ID 9360
Human Swiss-Prot Number Q13427

Alternative Names Peptidyl-prolyl cis-trans isomerase G (PPlase G)

(Peptidyl-prolyl isomerase G) (EC 5.2.1.8) (CASP10)

(Clk-associating RS-cyclophilin) (CARS-Cyp)


(CARS-cyclophilin) (SR-cyclophilin) (SR-cyp) (SRcyp)

(Cyclophilin G) (Rotamase G)

Background







Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950