

PP13G rabbit pAb

Cat No.: ES14015

For research use only

Overview

Product Name PP13G rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human PP13G AA

range: 11-61

Specificity This antibody detects endogenous levels of PP13G

at Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name PP13G
Gene Name PPP1R3G
Cellular localization cytoplasm,

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 648791 Human Swiss-Prot Number B7ZBB8


Alternative Names

+86-27-59760950

Background







Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night

