

## **POTEC** rabbit pAb

Cat No.: ES14018

For research use only

## Overview

Product Name POTEC rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human POTEC AA

range: 265-315

**Specificity** This antibody detects endogenous levels of POTEC

at Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name POTEC

Gene Name POTEC A26B2 POTE18

**Cellular localization** 

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 388468 Human Swiss-Prot Number B2RU33


**Alternative Names** 

+86-27-59760950

**Background** 







Western blot analysis of lysates from Hela cells, primary antibody was diluted at 1:1000, 4° over night

