

PDZD4 rabbit pAb

Cat No.: ES14171

For research use only

Overview

Product Name PDZD4 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1:500-2000

Immunogen Synthesized peptide derived from human PDZD4 AA

range: 4-54

Specificity This antibody detects endogenous levels of PDZD4

at Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name PDZD4

Gene Name PDZD4 KIAA1444 PDZK4 PDZRN4L

Cellular localization Cytoplasm, cell cortex . Mainly localized under the

plasma membrane.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 85kD
Human Gene ID 57595
Human Swiss-Prot Number Q76G19

Alternative Names PDZ domain-containing protein 4 (PDZ

domain-containing RING finger protein 4-like


protein)

Background

+86-27-59760950







Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

