

## **MYCT** rabbit pAb

## Cat No.:ES14637

For research use only

## Overview

MYCT rabbit pAb
Rabbit
WB;IHC
Human; Mouse;Rat
WB 1:500-2000;IHC-p 1:50-300
Synthesized peptide derived from human MYCT AA range: 352-402
This antibody detects endogenous levels of MYCT at Human/Mouse/Rat
Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Store at -20°C. Avoid repeated freeze-thaw cycles.
МҮСТ
SLC2A13
Cell membrane ; Multi-pass membrane protein .
The antibody was affinity-purified from rabbit
antiserum by affinity-chromatography using
epitope-specific immunogen.
Polyclonal
1 mg/ml
114134
Q96QE2




ELKbio@ELKbiotech.com

www.elkbiotech.com

www.eikbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night

Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C