

MXRA8 rabbit pAb


Cat No.:ES14640

For research use only


Overview

| | |
|--------------------------|---|
| Product Name | MXRA8 rabbit pAb |
| Host species | Rabbit |
| Applications | WB;IHC |
| Species Cross-Reactivity | Human; Mouse;Rat |
| Recommended dilutions | WB 1:500-2000;IHC-p 1:50-300 |
| Immunogen | Synthesized peptide derived from human MXRA8 AA range: 42-92 |
| Specificity | This antibody detects endogenous levels of MXRA8 at Human/Mouse/Rat |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Protein Name | MXRA8 |
| Gene Name | MXRA8 |
| Cellular localization | Cell membrane ; Single-pass type I membrane protein . Cell junction, tight junction . Cytoplasm . Cell projection, cilium membrane . Nucleus . Primarily localizes to the cell membrane (PubMed:23386276). Detected in the cilium of primary chondrocytes (By s |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 54587 |
| Human Swiss-Prot Number | Q9BRK3 |
| Alternative Names | |
| Background | |





Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).

