

Lyn (phospho Tyr508) rabbit pAb

Cat No.:ES1474

For research use only

Overview

Product Name	Lyn (phospho Tyr508) rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human Lyn around the
	phosphorylation site of Tyr507. AA range:463-512
Specificity	Phospho-Lyn (Y508) Polyclonal Antibody detects
opeometry	endogenous levels of Lyn protein only when
	phosphorylated at Y508.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Tyrosine-protein kinase Lyn
Gene Name	LYN
Cellular localization	Cell membrane. Nucleus. Cytoplasm. Cytoplasm,
	perinuclear region. Golgi apparatus. Membrane ;
	Lipid-anchor . Accumulates in the nucleus by
	inhibition of CRM1-mediated nuclear export.
	Nuclear accumulation is increased by inhibition of
	its kinase activity. The trafficking from the Golgi
	apparatus to the plasma membrane occurs in a
	kinase domain-dependent but kinase activity
	independent manner and is mediated by exocytic
	vesicular transport. Detected on plasma membrane
	lipid rafts.
Purification	The antibody was affinity-purified from rabbit
Fulfication	antiserum by affinity-chromatography using
	, , , , , , , , , , , , , , , , , , , ,
Clonality	epitope-specific immunogen. Polyclonal
_	roiycionai
	El Khie @El Khietech com



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Concentration	1 mg/ml
Observed band	53,56kD
Human Gene ID	4067
Human Swiss-Prot Number	P07948
Alternative Names	LYN; JTK8; Tyrosine-protein kinase Lyn;
	Lck/Yes-related novel protein tyrosine kinase;
	V-yes-1 Yamaguchi sarcoma viral related oncogene
	homolog; p53Lyn; p56Lyn
Background	This gene encodes a tyrosine protein kinase, which
	maybe involved in the regulation of mast cell
	degranulation, and erythroid differentiation.
	Alternatively spliced transcript variants encoding
	different isoforms have been found for this gene.
	[provided by RefSeq, Jul 2011],



Western Blot analysis of various cells using Phospho-Lyn (Y508) Polyclonal Antibody diluted at 1:1000

Western Blot analysis of 293 cells using Phospho-Lyn (Y508) Polyclonal Antibody diluted at 1:1000



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using Lyn (Phospho-Tyr507) Antibody

Immunohistochemistry analysis of paraffin-embedded human brain, using Lyn (Phospho-Tyr507) Antibody. The picture on the right is blocked with the phospho peptide.





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C