

MFSD7 rabbit pAb

Cat No.:ES14952

For research use only

Overview


Product Name	MFSD7 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human MFSD7 AA range: 259-309
Specificity	This antibody detects endogenous levels of MFSD7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	MFSD7
Gene Name	MFSD7 LP2561 UNQ385/PRO717
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kD
Human Gene ID	84179
Human Swiss-Prot Number	Q6UXD7
Alternative Names	Major facilitator superfamily domain-containing protein 7 (Myosin light polypeptide 5 regulatory protein) (MYL5)

Background





ELK Biotechnology



Western blot analysis of lysates from VEC cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.