



LRRT2 rabbit pAb

Cat No.:ES15095

For research use only


Overview

Product Name	LRRT2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human LRRT2 AA range: 372-422
Specificity	This antibody detects endogenous levels of LRRT2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	LRRT2
Gene Name	LRRTM2 KIAA0416 LRRN2
Cellular localization	Cell membrane; Single-pass type I membrane protein. Cell junction, synapse, postsynaptic cell membrane ; Single-pass type I membrane protein . Localized to excitatory synapses. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	26045
Human Swiss-Prot Number	O43300
Alternative Names	
Background	





ELK Biotechnology



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C