

LRIQ4 rabbit pAb

Cat No.:ES15107

For research use only


Overview

Product Name	LRIQ4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human LRIQ4 AA range: 480-530
Specificity	This antibody detects endogenous levels of LRIQ4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	LRIQ4
Gene Name	LRRIQ4 LRRC64
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kD
Human Gene ID	344657
Human Swiss-Prot Number	A6NIV6
Alternative Names	Leucine-rich repeat and IQ domain-containing protein 4 (Leucine-rich repeat-containing protein 64)
Background	





ELK Biotechnology



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.