



LLPH rabbit pAb

Cat No.:ES15144

For research use only


Overview

Product Name	LLPH rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human LLPH AA range: 18-68
Specificity	This antibody detects endogenous levels of LLPH at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	LLPH
Gene Name	LLPH C12orf31
Cellular localization	Nucleus, nucleolus . Chromosome . Cell-permeable protein. 22 hours after injection in the hippocampal area CA1, internalized by most cells at the injection site (By similarity). Localizes at the chromosome periphery during mitosis (PubMed:20813266). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84298
Human Swiss-Prot Number	Q9BRT6
Alternative Names	
Background	





ELK Biotechnology



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.