



LETM1 rabbit pAb

Cat No.:ES15167

For research use only

Overview

Product Name	LETM1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human LETM1 AA range: 574-624
Specificity	This antibody detects endogenous levels of LETM1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	LETM1
Gene Name	LETM1
Cellular localization	Mitochondrion inner membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3954
Human Swiss-Prot Number	O95202
Alternative Names	
Background	This gene encodes a protein that is localized to the inner mitochondrial membrane. The protein functions to maintain the mitochondrial tubular shapes and is required for normal mitochondrial morphology and cellular viability. Mutations in this gene cause Wolf-Hirschhorn syndrome, a complex malformation syndrome caused by the deletion of parts of the distal short arm of chromosome 4.





Related pseudogenes have been identified on chromosomes 8, 15 and 19. [provided by RefSeq, Oct 2009],

Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night

