

LACB2 rabbit pAb

Cat No.: ES15219

For research use only

Overview

Product Name LACB2 rabbit pAb

Host species Rabbit

Applications WB;ELISA;IHC

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000 **Immunogen** Synthesized peptide derived from human LACB2 AA

range: 19-69

Specificity This antibody detects endogenous levels of LACB2 at

Human/Mouse/Rat

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name LACB2

Gene Name LACTB2 CGI-83

Cellular localization Mitochondrion matrix .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 51110 Human Swiss-Prot Number Q53H82


Alternative Names

+86-27-59760950


Background







Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

