



# ITM2C rabbit pAb


Cat No.:ES15409

For research use only


## Overview

<b>Product Name</b>	ITM2C rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB;IHC
<b>Species Cross-Reactivity</b>	Human; Mouse;Rat
<b>Recommended dilutions</b>	WB 1:500-2000;IHC-p 1:50-300
<b>Immunogen</b>	Synthesized peptide derived from human ITM2C AA range: 1-51
<b>Specificity</b>	This antibody detects endogenous levels of ITM2C at Human/Mouse/Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	ITM2C
<b>Gene Name</b>	ITM2C BRI3 hucep-14 NPD018 PSEC0047
<b>Cellular localization</b>	Lysosome membrane ; Single-pass type II membrane protein . Cell membrane ; Single-pass type II membrane protein .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	81618
<b>Human Swiss-Prot Number</b>	Q9NQX7
<b>Alternative Names</b>	
<b>Background</b>	





Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

