

ABHD6 rabbit pAb

Cat No.: ES1573

For research use only

Overview

Product Name ABHD6 rabbit pAb

Host species Rabbit

Applications WB;IF;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions Western Blot: 1/500 - 1/2000. Immunofluorescence:

1/200 - 1/1000. ELISA: 1/20000. Not yet tested in

other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human ABHD6. AA

range:231-280

Specificity ABHD6 Polyclonal Antibody detects endogenous

levels of ABHD6 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Monoacylglycerol lipase ABHD6

Gene Name ABHD6

Cellular localization Late endosome membrane; Single-pass type II

membrane protein . Lysosome membrane ; Single-pass type II membrane protein .

Mitochondrion membrane; Single-pass type II

membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 38kD
Human Gene ID 57406
Human Swiss-Prot Number Q9BV23

Alternative Names ABHD6; Monoacylglycerol lipase ABHD6;

2-arachidonoylglycerol hydrolase; Abhydrolase

domain-containing protein 6




+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com




Background


catalytic activity:Hydrolyzes glycerol monoesters of long-chain fatty acids.,function:Has 2-arachidonoylglycerol hydrolase activity (By similarity). May be a regulator of endocannabinoid signaling pathways.,similarity:Belongs to the AB hydrolase superfamily.,




19-

Western Blot analysis of various cells using ABHD6 Polyclonal Antibody







Western blot analysis of lysates from COLO cells, using ABHD6 Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com