



GRAM2 rabbit pAb

Cat No.:ES15916

For research use only


Overview

Product Name	GRAM2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human GRAM2 AA range: 150-200
Specificity	This antibody detects endogenous levels of GRAM2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	GRAM2
Gene Name	GRAMD2
Cellular localization	Endoplasmic reticulum membrane ; Single-pass membrane protein . Cell membrane ; Peripheral membrane protein . Localizes to endoplasmic reticulum-plasma membrane contact sites (EPCS). Anchored at the ER, PM binding is mediated via GRAM domain and phosphati
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	38kD
Human Gene ID	196996
Human Swiss-Prot Number	Q8IU3
Alternative Names	GRAM domain-containing protein 2
Background	





ELK Biotechnology



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C