



# GNT2C rabbit pAb

Cat No.:ES16134

For research use only


## Overview

|                                 |   |
|---------------------------------|---|
| <b>Product Name</b>             | GNT2C rabbit pAb  |
| <b>Host species</b>             | Rabbit  |
| <b>Applications</b>             | WB  |
| <b>Species Cross-Reactivity</b> | Human;Rat;Mouse;  |
| <b>Recommended dilutions</b>    | WB 1: 500-2000  |
| <b>Immunogen</b>                | Synthesized peptide derived from human GNT2C AA range: 251-301  |
| <b>Specificity</b>              | This antibody detects endogenous levels of GNT2C at Human   |
| <b>Formulation</b>              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Storage</b>                  | Store at -20°C. Avoid repeated freeze-thaw cycles.  |
| <b>Protein Name</b>             | GNT2C   |
| <b>Gene Name</b>                | GCNT2 GCNT5 II NACGT1   |
| <b>Cellular localization</b>    |   |
| <b>Purification</b>             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Clonality</b>                | Polyclonal  |
| <b>Concentration</b>            | 1 mg/ml   |
| <b>Observed band</b>            |   |
| <b>Human Gene ID</b>            | 2651  |
| <b>Human Swiss-Prot Number</b>  | Q8NFS9  |
| <b>Alternative Names</b>        |   |
| <b>Background</b>               |   |





**ELK Biotechnology**



Western blot analysis of lysates from HuvEc cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

[ELKbio@ELKbiotech.com](mailto:ELKbio@ELKbiotech.com)

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.