

GNT2A rabbit pAb

Cat No.: ES16135

For research use only

Overview

Product Name GNT2A rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human;Rat;Mouse; Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human GNT2A AA

range: 191-241

Specificity This antibody detects endogenous levels of GNT2A

at Human

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name GNT2A

Gene Name GCNT2 GCNT5 II NACGT1

Cellular localization Golgi apparatus membrane ; Single-pass type II

membrane protein.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 2651 Human Swiss-Prot Number Q8N0V5


Alternative Names

+86-27-59760950

Background







Western blot analysis of lysates from Hela cells, primary antibody was diluted at 1:1000, 4° over night

