

## AK5 rabbit pAb

## Cat No.:ES1619

For research use only

## Overview

Product Name	AK5 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse
<b>Recommended dilutions</b>	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/5000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human AK5. AA range:391-440
Specificity	AK5 Polyclonal Antibody detects endogenous levels
	of AK5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Adenylate kinase isoenzyme 5
Gene Name	АК5
<b>Cellular localization</b>	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	25kD
Human Gene ID	26289
Human Swiss-Prot Number	Q9Y6K8
Alternative Names	AK5; Adenylate kinase isoenzyme 5; AK 5; ATP-AMP
	transphosphorylase 5
Background	adenylate kinase 5(AK5) Homo sapiens This
	gene encodes a member of the adenylate kinase
	family, which is involved in regulating the adenine
	nucleotide composition within a cell by catalyzing
	the reversible transfer of phosphate groups among
	adenine nucleotides. This member is related to the



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



UMP/CMP kinase of several species. It is located in the cytosol and expressed exclusively in brain. Alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using AK5 Polyclonal Antibody

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using AK5 Antibody. The picture on the right is blocked with the synthesized peptide.



Jurka

(kD)

117-85-

48-

34-

26-

19-



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C