

## FOXD2 rabbit pAb

Cat No.: ES16347

For research use only

## Overview

Product Name FOXD2 rabbit pAb

Host species Rabbit Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1: 500-2000

Immunogen Synthesized peptide derived from human FOXD2 AA

range: 177-227

**Specificity** This antibody detects endogenous levels of FOXD2

at Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name FOXD2

Gene Name FOXD2 FKHL17 FREAC9

Cellular localization Nucleus.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 2306 Human Swiss-Prot Number 060548

**Alternative Names** 


**Background** This gene belongs to the forkhead family of

transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined. [provided by

RefSeq, Jul 2008],







Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night

