



FDXA1 rabbit pAb

Cat No.:ES16492

For research use only


Overview

Product Name	FDXA1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human FDXA1 AA range: 358-408
Specificity	This antibody detects endogenous levels of FDXA1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	FDXA1
Gene Name	FDXACB1
Cellular localization	cytoplasm,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	91893
Human Swiss-Prot Number	Q9BRP7
Alternative Names	
Background	This gene encodes a protein which contains a ferredoxin-fold anticodon-binding domain which is contained in a subset of phenylalanyl tRNA synthetases. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2011],





ELK Biotechnology



Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.