



FAT3 rabbit pAb

Cat No.:ES16518

For research use only

Overview

Product Name	FAT3 rabbit pAb
Host species	Rabbit
Applications	IHC;IF
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	IHC-p 1: 50-200
Immunogen	Synthesized peptide derived from human FAT3 AA range: 896-946
Specificity	This antibody detects endogenous levels of FAT3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	FAT3
Gene Name	FAT3 CDHF15 KIAA1989
Cellular localization	Membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	120114
Human Swiss-Prot Number	Q8TDW7
Alternative Names	
Background	

