



FA72A rabbit pAb

Cat No.:ES16564

For research use only


Overview

Product Name	FA72A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human FA72A AA range: 85-135
Specificity	This antibody detects endogenous levels of FA72A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	FA72A
Gene Name	FAM72A UGENE
Cellular localization	Cytoplasm. Mitochondrion. A V5 epitope-tagged construct has been shown to localize to the nucleus (PubMed:18676834). 5-7% of total FAM72A is associated with mitochondria around the nucleus in HEK293 cells (PubMed:21317926). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	729533
Human Swiss-Prot Number	Q5TYM5
Alternative Names	
Background	





ELK Biotechnology



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.