

FA43A rabbit pAb


Cat No.:ES16577

For research use only

Overview

Product Name	FA43A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse
Recommended dilutions	WB 1: 500-2000
Immunogen	Synthesized peptide derived from human FA43A AA range: 230-280
Specificity	This antibody detects endogenous levels of FA43A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	FA43A
Gene Name	FAM43A PP7298
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	131583
Human Swiss-Prot Number	Q8N2R8
Alternative Names	
Background	





Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night

