

F16B1 rabbit pAb

Cat No.: ES16621

For research use only

Overview

Product Name F16B1 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse Recommended dilutions WB 1:500-2000

Immunogen Synthesized peptide derived from human F16B1 AA

range: 211-261

Specificity This antibody detects endogenous levels of F16B1 at

Human/Mouse

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name F16B1

Gene Name FAM160B1 KIAA1600

Cellular localization

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.


Clonality Polyclonal
Concentration 1 mg/ml
Observed band 85kD
Human Gene ID 57700
Human Swiss-Prot Number Q5W0V3

Alternative Names Protein FAM160B1

Background







Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

