



# ELOF1 Polyclona Antibody

Cat No.:ES16758

For research use only


## Overview

Product Name	ELOF1 Polyclona Antibody
Host species	Rabbit
Applications	IHC; ELISA
Species Cross-Reactivity	Human; Mouse; Rat
Recommended dilutions	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Immunogen	Synthesized peptide derived from human ELOF1 AA range: 30-110
Specificity	This antibody detects endogenous levels of human ELOF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ELOF1
Gene Name	ELOF1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84337
Human Swiss-Prot Number	P60002
Alternative Names	Transcription elongation factor 1 homolog
Background	function:Transcription elongation factor implicated in the maintenance of proper chromatin structure in actively transcribed regions.,similarity:Belongs to the ELOF1 family.,





**ELK Biotechnology**



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

[ELKbio@ELKbiotech.com](mailto:ELKbio@ELKbiotech.com)

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C