

eIF4B (phospho-Ser406) rabbit pAb

Cat No.:ES16769

For research use only

Overview

Product Name	eIF4B (phospho-Ser406) rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized phospho peptide around human eIF4B (Ser406)
Specificity	This antibody detects endogenous levels of Human Mouse Rat eIF4B (phospho-Ser406)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	eIF4B (Ser406)
Gene Name	EIF4B
Cellular localization	cytosol,polysome,eukaryotic translation initiation factor 4F complex,dendrite,neuronal cell body,
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	80kD
Human Gene ID	1975
Human Swiss-Prot Number	P23588
Alternative Names	Eukaryotic translation initiation factor 4B (eIF-4B)
Background	

