



DPYL3 rabbit pAb

Cat No.:ES16904

For research use only

Overview


Product Name	DPYL3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human; Mouse;Rat
Recommended dilutions	WB 1:500-2000
Immunogen	Synthesized peptide derived from human DPYL3 AA range: 240-290
Specificity	This antibody detects endogenous levels of DPYL3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	DPYL3
Gene Name	DPYSL3 CRMP4 DRP3 ULIP ULIP1
Cellular localization	Cytoplasm . Cell projection, growth cone . Colocalizes with synaptic vesicle protein 2 in the central region of the growth cone. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	65kD
Human Gene ID	1809
Human Swiss-Prot Number	Q14195
Alternative Names	Dihydropyrimidinase-related protein 3 (DRP-3) (Collapsin response mediator protein 4) (CRMP-4) (Unc-33-like phosphoprotein 1) (ULIP-1)

Background





ELK Biotechnology



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C