



# CRAS1 rabbit pAb

Cat No.:ES17234

For research use only

## Overview

<b>Product Name</b>	CRAS1 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	IHC;IF
<b>Species Cross-Reactivity</b>	Human;Rat;Mouse;
<b>Recommended dilutions</b>	IHC-p 1: 50-200
<b>Immunogen</b>	Synthesized peptide derived from human CRAS1 AA range: 54-104
<b>Specificity</b>	This antibody detects endogenous levels of CRAS1 at Human
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	CRAS1
<b>Gene Name</b>	CRYM-AS1 NCRNA00169
<b>Cellular localization</b>	Membrane ; Single-pass membrane protein .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	
<b>Human Swiss-Prot Number</b>	A6NIL9
<b>Alternative Names</b>	
<b>Background</b>	

