



CDC20 (phospho-Ser51) rabbit pAb

Cat No.:ES17568

For research use only

Overview

Product Name	CDC20 (phospho-Ser51) rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized phospho peptide around human CDC20 (Ser51)
Specificity	This antibody detects endogenous levels of Human Mouse CDC20 (phospho-Ser51)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	CDC20 (Ser51)
Gene Name	CDC20
Cellular localization	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle pole .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50kD
Human Gene ID	991
Human Swiss-Prot Number	Q12834
Alternative Names	Cell division cycle protein 20 homolog (p55CDC)
Background	CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation. [provided by RefSeq, Jul 2008],

